
Deutsch	4	Nederlands	16
CERABLATE® easy/easy TC Steuerbarer Ablationskatheter		CERABLATE® easy/easy TC Bestuurbare ablatiekatheter	
Gebrauchsanweisung		Gebruiksaanwijzing	
English	5	Polski	17
CERABLATE® easy/easy TC Steerable Ablation Catheter		CERABLATE® easy/easy TC Sterowalny cewnik ablacyjny	
Instructions for Use		Instrukcja użycia	
Français	6	Português	18
CERABLATE® easy/easy TC Cathéter d'ablation orientable		CERABLATE® easy/easy TC Cateter de ablação controlável	
Instructions d'utilisation		Instruções de operação	
Italiano	8	Svenska	20
CERABLATE® easy/easy TC Catetere per ablazione manovrabile		CERABLATE® easy/easy TC Styrbar ablationskateter	
Istruzioni per l'uso		Bruksanvisning	
Español	9	Slovensky	21
CERABLATE® easy/easy TC Catéter de ablación dirigible		CERABLATE® easy/easy TC Ovládatelný ablačný katéter	
Instrucciones de uso		Návod na použitie	
Dansk	10	Slovenščina	22
CERABLATE® easy/easy TC Styrbar ablationskateter		CERABLATE® easy/easy TC Krmiljeni ablačijski kateter	
Brugsanvisning		Navodila za uporabo	
Suomi	12	Čeština	23
CERABLATE® easy/easy TC Ohjattava ablaatio-katetri		CERABLATE® easy/easy TC Ovládatelný ablační katetr	
Käyttöohje		Návod k použití	
Ελληνικά	13	Türkçe	25
CERABLATE® easy/easy TC Ελεγχόμενος καθετήρας εκτομής		CERABLATE® easy/easy TC Yönlendirilebilir Ablasyon Kateteri	
Οδηγίες χρήσης		Kullanım Kılavuzu	
Latviski	14	Magyar	26
CERABLATE® easy/easy TC Vadāms ablācijas katetrs		CERABLATE® easy/easy TC vezérelhető ablációs katéter	
Lietošanas instrukcija		Használati útmutató	

	de • Herstellungsdatum en • Date of manufacture fr • Date de fabrication it • Data di costruzione es • Fecha de fabricación da • Produktionsdato fi • Valmistuspäivämäärä el • Ημερομηνία κατασκευής lv • Izlaiduma datums nl • Fabricagedatum pl • Data produkcji pt • Data de fabrico sv • Tillverkningsdatum sk • Dátum výroby sl • Datum proizvodnje es • Datum výroby tr • Üretim tarihi hu • Gyártási dátum		de • Sterilisation mit Ethylenoxid en • Sterilisation with ethylene oxide fr • stérilisation à l'oxyde d'éthylène it • Sterilizzazione con ossido di etilene es • Esterilizado con óxido de etileno da • Steriliseret med ethylenoxidgas fi • Steriloitu etyleenioksidikaasulla el • Αποστείρωση με αέριο αιθυλενοξείδιο lv • Sterilizēts ar etilēna oksīda gāzi nl • Gesteriliseerd met ethylenoxidegas pl • Wysterylizowane tlenkiem etylenu pt • Esterilizado com óxido de etileno sv • Steriliserad med etylenoxidgas sk • Sterilizované etylénoxidom sl • Sterilizirano z etilenoksidom es • Sterilizováno ethylenoxidovým plynem tr • Etilen oksit gazı ile sterilize edildi hu • Etilénoxid gázzal sterilizálva		de • Trocken aufbewahren en • Keep dry fr • Conserver au sec it • Mantenere secco es • Mantener seco da • Opbevares tørt fi • Suojaa märältä el • Διατηρείται στεγνό lv • Sargāt no mitruma nl • Tegen vocht beschermen pl • Chronić przed wilgocią pt • Manter seco sv • Förvaras tørt sk • Chránit' pred vlhkom sl • Varujte pred vlago es • Cbrañte pñed mokrem tr • Rutubetten koruyunuz hu • Szárazon tartandó!
	de • Hersteller en • Manufacturer fr • Fabricant it • Fabbricante es • Fabricante da • Producent fi • Valmistaja el • Κατασκευαστής lv • Ražotājs nl • Fabrikant pl • Producent pt • Fabricante sv • Tillverkare sk • Výrobca sl • Proizvajalec es • Výrobce tr • Üretici hu • gyártó		de • Gebrauchsanweisung beachten en • Consult Instructions for Use fr • Consulter les instructions d'utilisation it • Consultare le istruzioni per l'uso es • Consulte las instrucciones de uso da • Følg brugsanvisningen fi • Noudata käyttöohjetta el • Προσέχετε τις οδηγίες χρήσης lv • Ievērot lietošanas instrukciju nl • Gebruiksaanwijzing raadplegen pl • Przestrzegać instrukcji użytkowania pt • Observar as instruções de operação sv • Beakta bruksanvisningen sk • Rešpektujte návod na použitie sl • Upoštečajte navodila za uporabo es • Respetujte návod k použití tr • Kullanım talimatlarını izleyin hu • tartsa be a használati utasítást		de • Bogendurchmesser en • Arc diameter fr • Diamètre d'arc it • Diametro arco es • Diámetro de la curva da • Buediameter fi • Kaaren halkaisija el • Διάμετρος τόξου lv • Loka diametrs nl • Boogdiameter pl • Średnica łuku pt • Diâmetro do arco sv • Krökningsdiameter sk • Priemer oblúka sl • Premer loka es • Príměr oblouku tr • Kavis çapı hu • ív átmérője
	de • Nur zum Einmalgebrauch en • For one-time usage only fr • A usage unique it • Monouso es • De usar sólo una vez da • Kun til engangsbrug fi • Kertäkäyttötuote el • Για μία μόνο χρήση lv • Tikai vienreizējai lietošanai nl • Uitsluitend bestemd voor eenmalig gebruik pl • Tylko do jednorazowego użytku pt • Utilizar apenas uma vez sv • Endast för engångsbruk sk • Len na jedno použitie sl • Samo za enkratno uporabo es • Pouze k jednorázovému použití tr • Sadece tek kullanımlık hu • Csak egyszeri használatra		de • Bestellnummer en • Order number fr • Référence de commande it • Numero d'ordine es • Número de referencia da • Bestillingsnummer fi • Tilausnumero el • Αριθμός παραγγελίας lv • Pasūtījuma numurs nl • Bestelnummer pl • Numer zamówienia pt • Número de encomenda sv • Ordernummer sk • Číslo objednávky sl • Številka za naročanje es • Číslo objednávky tr • Sipariş numarası hu • Megrendelési szám		de • Kurve en • Curvature fr • Courbe it • Curva es • Curva de • Kurve fi • Käyrä el • Καμπύλη lv • Līkne nl • Curve pl • Krzywa pt • Curva sv • Kurva sk • Krivka sl • Krivulja es • Křivka tr • Eğri hu • görbület
	de • Chargen-Bezeichnung en • Batch name fr • Désignation des charges it • Lotto, designazione es • Descripción de la carga da • Lotbetegnelse fi • Valmistuserän kuvaus el • Περιγραφή παρτίδας lv • Partijas apzīmējums nl • Lotaanduiding pl • Nazwa partii pt • Designação do lote sv • Partibeteckning sk • Označenie šarže sl • Oznaka šarže es • Číslo šarže tr • Parti adı hu • Gyártási tétele szám		de • Inhalt beschädigter Packung nicht verwenden en • Do not use if package damaged fr • Ne pas utiliser si l'emballage est endommagé it • Non utilizzare se la confezione è danneggiata es • No usar si el paquete está dañado da • Indhold fra beskadiget emballage må ikke anvendes fi • Älä käyttää vaurioituneen pakkauksen sisältöä el • Μη χρησιμοποιείτε το περιεχόμενο φθαρμένου συσκευασιστών lv • Neizmantot bojāta iepakojuma saturu nl • Inhoud van een beschadigde verpakking niet gebruiken pl • Nie używać w przypadku uszkodzonego opakowania pt • Não usar o conteúdo de embalagem danificada sv • Använd inte innehålllet vid skadad förpackning sk • Obsah poškodených obalov nepoužívajte sl • Ne uporabite vsebine v poškodovani embalaži es • Obsah poškozeného obalu nepoužívejte tr • Hasarlı paket içeriğini kullanmayın hu • ha a csomagolás megrongálódott, a tartalmát ne használja fel		de • Benutzbare Länge en • Usable length fr • Longueur utile it • Lunghezza utile es • Longitud utilizable da • Anvendelseslængde fi • Käytettävä pituus el • Χρήσιμο μήκος lv • Lietojamais garums nl • Bruikbare lengte pl • Długość użytkowa pt • Comprimento útil sv • Effektiv längd sk • Použitelná dĺžka sl • Uporabna dolžina es • Pouzitelná délka tr • Kullanılabilir uzunluk hu • használható hossz
	de • Verwendbar bis en • To be used before fr • Utiliser avant it • Utilizzabile fino a es • Utilizar antes da • Holdbar til fi • Viimeinen käyttöpäivä el • Χρήση έως lv • Lietojams līdz nl • Te gebruiken tot pl • Przydatne do użycia do pt • Utilizar até sv • Används före sk • Použiteľné do sl • Uporabno do es • Spotřebujte do tr • Son kullanma tarihi hu • Minőségét megőrzi		de • Nicht erneut sterilisieren en • Do not sterilise anew fr • Ne pas stériliser à nouveau it • Non sterilizzare di nuovo es • No esterilizar de nuevo da • Må ikke gensteriliseres fi • Ei saa steriloida uudestaan el • Μην επανααστείρωσετε lv • Nesterilizēt atkārtoti nl • Niet opnieuw steriliseren pl • Nie sterylizować ponownie pt • Não esterilizar novamente sv • Får ej resteriliseras sk • Opakovane nesterilizujte sl • Ne sterilizirajte znova es • Opakovaně nesterilizujte tr • Tekrar sterilize etmeyin hu • Ne sterilizálja újból		de • Elektrodenanzahl en • Number of electrodes fr • Nombre d'électrodes it • Numero di elettrodi es • Número de electrodos da • Elektrodeantal fi • Elektrodien lukumäärä el • Αριθμός ηλεκτροδίων lv • Elektrodu skaits nl • Aantal elektroden pl • Liczba elektrod pt • Quantidade de eletrodos sv • Elektrodeantal sk • Počet elektród sl • Število elektrod es • Počet elektrod tr • Elektrot sayısı hu • elektródák száma
	de • Temperaturbegrenzung en • Temperature limitation fr • Limites de température it • Limiti di temperatura es • Limite de temperatura da • Temperaturbegrænsning fi • Lämpötilan rajoitus el • Οριοθέτηση θερμοκρασίας lv • Temperatūras ierobežojums nl • Temperatuurgrens pl • Ograniczenie temperatury pt • Limite de temperatura sv • Temperaturbegränsning sk • Teplotné ohraničenie sl • Omejitve temperature es • Teplotní omezení tr • Sıcaklık limiti hu • Hőmérséklethatárok		de • Vor Sonnenlicht geschützt aufbewahren en • Keep away from Sunlight fr • Ne pas exposer aux rayons solaires it • Evitare l'esposizione alla luce del sole es • Evitar la exposición a los rayos solares da • Må ikke udsættes for sollys fi • Suojaa auringonvalolta el • Προστατέψτε από την ηλιακή ακτινοβολία lv • Sargāt no saules stariem nl • Tegen zonlicht beschermen pl • Chronić przed nasłonecznieniem pt • Proteger da luz solar sv • Skyddas mot direkt solljus sk • Chránit' pred slnečným svetlom sl • Varujte pred sončno svetlobo es • Cbrañte pñed slnečným svĕtlem tr • Güneş ışınlarından koruyunuz hu • Védjük a napfénytől!		de • Elektrodenabstand en • Electrode spacing fr • Écart des électrodes it • Distanza elettrodi es • Distancia entre electrodos da • Elektrodeafstand fi • Elektrodien etäisyys el • Απόσταση ηλεκτροδίων lv • Elektrodu atālums nl • Elektrode-afstand pl • Odstęp między elektrodami pt • Distância dos eletrodos sv • Elektrodavstånd sk • Vzdialenosť elektród sl • Razdalja med elektrodami es • Rozestup elektrod tr • Elektrot mesafesi hu • elektródák távolsága

	de • Elektrodenmaterial en • Electrode material fr • Matériau de l'électrode it • Materiale elettrodo es • Material del electrodo da • Elektrodemateriale fi • Elektrodin materiaali el • Υλικό ηλεκτρόδιου lv • Elektroda materiāls nl • Elektrodemateriaal pl • Materiał elektrody pt • Material do eletrodo sv • Elektrodmaterial sk • Materiál elektródy sl • Material elektrod cs • Materiál elektrody tr • Elektrot malzemesi hu • elektróda anyaga		de • Elektrodenlänge en • Electrode length fr • Longueur de l'électrode it • Lunghezza elettrodo es • Longitud del electrodo da • Elektrodelængde fi • Elektrodin pituus el • Μήκος ηλεκτρόδιου lv • Elektroda garums nl • Elektrodelengte pl • Długość elektrody pt • Comprimento do eletrodo sv • Elektrodlängd sk • Dĺžka elektródy sl • Dolžina elektrod cs • Délka elektrody tr • Elektrot uzunluğu hu • elektróda hossza		de • Temperatursensor en • Temperature sensor fr • Capteur de température it • Sensore termico es • Sensor de temperatura da • Temperatursensor fi • Lämpötila-anturi el • Αισθητήρας θερμοκρασίας lv • Temperatūras devējs nl • Temperatuursensor pl • Czujnik temperatury pt • Sensor de temperatura sv • Temperatursensor sk • Snímač teploty sl • Temperaturno tipalo cs • Snímač teploty tr • Sıcaklık sensörü hu • hőmérséklet-érzékelő
	de • Durchmesser en • Diameter fr • Diamètre it • Diametro es • Diámetro da • Diameter fi • Halkaisija el • Διάμετρος lv • Diametrs nl • Diameter pl • Średnica pt • Diâmetro sv • Diameter sk • Priemer sl • Premer cs • Průměr tr • Çap hu • átmérő		de • Anschlussskabel en • Connector cable fr • Câble de connexion it • Cavo di connessione cavo es • Cable de conexión, cable da • Forbindelseskabel fi • Liitintäkkäpeli kaapeli el • Καλώδιο σύνδεσης lv • Savienotājkabelis nl • Aansluitkabel kabel pl • Przyłącze Kabel pt • Cabode conexão sv • Anslutningskabel sk • Pripojný kábel sl • Povezovalni kabel cs • Pripojný kabel tr • Bağlantı kablosu kablo hu • csatlakozókábel kábel		de • Inhalt en • Contents fr • Sommaire it • Contenuto es • Contenido da • Indhold fi • Sisältö el • Περιεχόμενα lv • Saturs nl • Inhoud pl • Treść pt • Conteúdo sv • Innehåll sk • Obsah sl • Vsebina cs • Obsah tr • İçindekiler hu • Tartalom
TC	de • Thermocouple en • Thermocouple fr • Thermocouple it • Termocoppia es • Termopar da • Thermocouple fi • Thermocouple el • Θερμοζεύγος lv • Termopāris nl • Thermocouple pl • Termoogniwo pt • Termopar sv • Thermocouple sk • Thermocouple sl • Termočlen cs • Thermocouple tr • Isıl çift hu • Thermocouple	TM	de • Thermistor en • Thermistor fr • Thermistance it • Termistore es • Termistor da • Termistor fi • Termistori el • Θερμίστορ lv • Termistors nl • Thermistor pl • Termistor pt • Termistor sv • Termistor sk • Termistor sl • Termistor cs • Termistor tr • Termistör hu • Termisztor		

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Steerable Ablation Catheter

General notes

Caution: Please read these instructions carefully before using the product.

1. Inspection before use

The product, packaging and any accessories must be inspected for damage and malfunction before use. The product must not be used if any damage or malfunction is encountered. If the instructions for use are no longer available you can request replacement instructions from OSYPKA AG. The instructions for use may only be disposed of after all products in the product unit have been used.

2. Use and handling

Caution: Use of the product in sterile conditions.

The product may only be used in a medical treatment facility that has been specifically set up for the appropriate application and by trained staff. The product may only be used if its safe application can be guaranteed. It is the doctor's responsibility to select the medically appropriate procedure and method. The operating instructions are intended purely as general information on handling the product. The product must not be modified in any way.

Warnings / general information / precautionary measures must always be observed even if a differing procedure is chosen for medical reasons! Failure to observe these can result in complications and malfunction!

3. Sterilisation and storage

The products are sterilised with ethylene oxide. Provided the sterilised product is stored correctly, it may be used up until expiry date (use before date) indicated on the package. Store the steril product in its original packaging in a cool, dry place away from light at between 10°C and 25°C.

Caution: Not autoclavable. Single use only. Do not re-sterilise used product. Re-use brings the risk of infection and the possibility of malfunction!

4. Disposal

The product may be contaminated after use. For this reason please observe your hospital's regulations regarding disposal. Any pointed or sharp objects such as needles, hypodermic needles or scalpels that were part of the product must be disposed of safely and separately in order to prevent injury and infection.

5. Liability limitation

Medical products from OSYPKA AG are manufactured from high-quality materials using managed and proven production processes. Quality is continuously monitored during manufacture and prior to delivery. Should you nevertheless encounter a functionally impaired or inoperative product, please return it once the defect occurs indicating the nature of the defect and impairment. We will then inspect the returned product immediately. In the event of a justified complaint we will provide a replacement free of charge.


Product damage that results from improper storage, handling, use for which the product is not intended, unauthorised modification to the product or prohibited re-use and re-sterilisation is excluded from warranty. The warranty will become void if the product is interfered with or modified by persons who have not been authorised to do so in writing by the manufacturer.

1. Product specification

CERABLATE® catheters were specially developed for intracardial radio-frequency ablation in connection with the ablation system HAT300 Smart or equivalent. The temperature at the catheter tip created by RF dissipation is continuously measured by a temperature sensor contained in the catheter tip. This enables accurate temperature control and precise temperature monitoring. Depending on the selected target value and the currently measured temperature, a microprocessor, integrated in the ablation devices, continuously controls RF dissipation. This prevents overheating or tissue carbonization. The position of thermo elements in the catheter tip provides excellent heat transfer and thus an optimum temperature profile during RF dissipation.

The arrangement of the individual poles (2-pin or 4-pin) enables precise recording of the intracardial signals during mapping and RF dissipation. The catheter tip ensures an adequate surface for effective ablation. CERABLATE® is a catheter with highly sensitive controllability, combined with excellent torsional control. This enables accurate catheter positioning. The inner pull wire, manufactured from a special material, ensures the catheter tip's extraordinary flexibility, which enables accurate adaptation of the distal catheter part to the irregular anatomy of the heart. The curve of the catheter tip can be infinitely adjusted by activating the manipulator at the handle of the CERABLATE® catheter up to at least 180°. This specially developed manipulator allows the variation of the curvature and maintaining its stability in every desired position. In models with the appropriate label or marking on the manipulator, the values of the manipulation forces can be adjusted.


Distal curvature control

Infinitely adjustable catheter curvature is achieved by moving the front part of the plastic handle (manipulator) forwards and backwards.

In models where the manipulation forces cannot be adjusted, the curvature of the tip is latched in each position. In catheters having adjustable manipulation forces, you have to proceed as follows for the adjustment:

Adjusting the manipulation forces

Hold the handle tight and turn the manipulator according to the label or marking. Turning in the direction of the open lock icon reduces the manipulation forces. Turning in the direction of the closed lock icon increases the manipulation forces and latches the curve of the tip. The manipulator should be rotated by no more than 1 turn in each direction from the factory setting!

Models:



CERABLATE® easy / easy TC

Standard RF catheter

Curve form: Classic

Curve form: Sigma

Temperature sensor:
Thermistor / thermoelement (thermocouple)

2. Indication

Radio-frequency catheter ablation with CERABLATE® catheters is used for:

- Supraventricular tachycardia
- Ventricular tachycardia
- Atrio- / ventricular tachycardia

3. Handling

The catheter should only be used by electrophysiologists who have been trained in intracardial radio-frequency catheter ablation and who are practiced in the technique.

The following notes are for informational purposes only and should be supplemented by the electrophysiologist on the basis of his/her clinical experience. National or international guidelines are to be considered.

- The catheters must be removed from their packaging under sterile conditions and used in a sterile work area.
- Insert the catheter into the vascular system via normal accesses for catheter systems.
- Connect the catheter to the radio-frequency ablation device (see operating instructions for the used device) or to the ECG recording device via suitable connection cables.
- Insert the catheter into the endocardial area using the X-ray and ECG system.
- Before pulling out the catheter from the patient, it must be checked with the help of X-ray that the tip of the catheter is not bent!
- For better positioning in the heart, you can bend the catheter tip with the use of the control mechanism and thus change the catheter curvature. Information on the control mechanism is contained in Chapter 1 „Product description“.
- You can start RF dissipation once the lead tip displays stable tissue contact with the intended ablation site.
- RF dissipation can be repeated at the same or other sites. Should the impedance increase for any reason, remove the catheter from the vessel. Clean the lead tip before repeating RF dissipation.

Cleaning:

Should you need to clean the catheters for any reason, be sure to protect the handle from penetrating fluids.

Caution: Liquid in the handle leads to an electrical short circuit and to the destruction of the catheter.

4. Patients

CERABLATE®: unlimited

5. Warnings

- Warning:** The catheter represents a direct, low-resistance power path to the heart. It is known that small electric currents are sufficient to make the heart flutter. Do not touch the connector with your bare hands and do not allow it to make contact with electrically conducting or damp surfaces.
- Never touch the catheter tip and the indifferent lead at the same time, especially not during output.
 - Keep away any static electricity from the catheter system.
 - When connecting the catheter to the ablation device, ensure compatibility of the connections with the ablation device. The catheter should only be used if the temperature display of the ablation device displays a realistic temperature.
 - Observe the respective operating instructions when connecting and using external devices.
 - A working external pacemaker and defibrillator should be available before and during the entire electrophysiological examination. During defibrillation, defibrillation discharges can result in current flows through the pacemaker/catheter circuit. To protect the patient, the catheter circuit must therefore be interrupted during defibrillation.

Caution: If devices operated with mains voltage are used near the patient, follow the instructions for electromedical devices to avoid creeping currents to the heart.

6. Precautionary measures

It is imperative to carefully check the catheter and its accessories in connection with the radio-frequency ablation system before using it.

Pay particular attention to whether the catheter shaft and the lead tip have been damaged in any way. Replace faulty products.

In catheters where the catheter tip can be bent, the working of the manipulator should be tested before beginning. (see 1. Product specification).

Caution: In models with adjustable manipulation forces, the manipulator should not be turned too far against the specifications given in the product description in the direction of the open lock, since, otherwise, the buckling protection may get shifted and thus become ineffective. Products affected in this manner must be replaced immediately.

All improper use before or during the application, such as excessive bending or kinking of the catheter, can lead to the destruction of the control mechanism.

The catheter should only be applied by specialists trained in the method of intracardial radio-frequency catheter ablation/electrophysiology and practiced in the technique.

- The function of active implantable medical products such as pacemakers and defibrillators can be impaired by radio-frequency catheter ablation.
- Radio-frequency catheter ablation causes electrical discharges in the catheter tip/tool tip, which may influence the function of external devices. Observe the operating instructions of the surrounding devices.
- The procedure must take place under X-ray conditions with appropriate X-ray equipment.
- The catheter's function can be influenced by magnetic fields.

Caution: Do not apply the catheter during MRI examinations.

- The application of external devices emitting electromagnetic radiation can lead to interference.
- The function of cardiological implants can be impaired by catheter/lead placement. Such systems could further impair the function of the ablation catheter. Take the patient's medical history into consideration.
- Unforeseen complications may occur during transseptal access.
- Observe the operating instructions of the used radio-frequency generator/ablation generator.
- Do not clean the handle with liquids. This could lead to a short circuit.
- If RF emission is output-controlled, this may lead to an increased complication rate. Thrombus formation and/or carbonization of the tissue in particular may increase. The possibility of the impairment of catheter function caused by overheating cannot be excluded in this case.

Caution: When using a defibrillator, disconnect the ablation catheter from the connection cable or remove the catheter from the endocardium. Do not defibrillate via the ablation catheter!

- The CERABLATE® catheter was specially developed for intracardial radio-frequency ablation in connection with the ablation system HAT300Smart.

Caution: The use of external devices may cause unforeseen risks. Clarify these with the device manufacturer.

When using the products with external generators, peak or effective voltage may not exceed 440 V_{ss} or 150 V_{eff} (at 600 Ω). RF impedance may not be less than 40 Ω.

- When connecting CERABLATE® catheters to other ablation systems, first check for compliance of PIN allocation. If compliance with the PIN allocation cannot be guaranteed, CERABLATE® may not be connected to other ablation systems.

Caution: The products should only be operated in temperature-controlled mode. Use the output limit and impedance disconnection functions.

7. Contraindication

- Patients after acute myocardial infarction
- Systemic infections
- Acute local infections at puncture site
- Patient has shown sensitivity to foreign objects in the past or has an extreme allergy
- Insufficient or unsuitable tissue, e.g. tissue damaged by radiation or suppuration
- Physiological or anatomical abnormalities which could lead to post-operative complications, such as bleeding diathesis, reduced resistance against infections, etc.
- Hemodynamic instability
- Patients unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation
- Patients with prosthetic valves

8. Possible complications / side effects

- Infection / sepsis
- Embolisms
- Coagulum formation, in particular with the use of catheters in power-regulated mode
- Arterial spasms (including coronary arteries)
- Perforations
- Pericardial effusion
- Pericardial tamponade
- Endocarditis
- AV-block
- Abnormalities of cardiac wall or heart valves
- Thrombosis of arteries or veins
- Obstruction of the left coronary artery
- Severe cardiac arrhythmias
- Ventricular fibrillation
- Stroke
- Death

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