CE 2797

Conformité Européenne (European Conformity), MDR 2017/745
Conformité Européenne (Europäische Konformität), MDR 2017/745
Conformité Européenne (Conformité Européenne), MDR 2017/745
Conformité Européenne (Conformità europea), MDR 2017/745
Conformité Européenne (Conformidad Europea), MDR 2017/745
Conformité Européenne (europaisk overensstemmelse) MDR 2017/745
Conformité Européenne (eurooppalainen vaatimustenmukaisuus), MDR 2017/745
Conformité Européenne (Ευρωπαϊκή Συμμόρφωση), Κανονισμός 2017/745 για τα Ιατροτεχνολογικά Προϊόντα (MDR)
Conformité Européenne (Eiropas atbilstība) MDR 2017/745
Conformité Européenne (Europese conformiteit), MDR 2017/745
Conformité Européenne (Zgodność europejska), MDR 2017/745
Conformité Européenne (Conformidade Europeia), MDR 2017/745
Conformité Européenne (Europeisk överensstämmelse), MDR 2017/745
Conformité Européenne (Európske orgány pre posudzovanie zhody), MDR 2017/745
Conformité Européenne (Evropska skladnost), MDR 2017/745
Conformité Européenne (Evropská shoda), MDR 2017/745
Conformité Européenne (Avrupa Uygunluğu), MDR 2017/745
Conformité Européenne (Európai megfelelés), MDR 2017/745

UK
CA 0086

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
CERABLATE® easy/easy TC

Steerable Ablation Catheter

Caution: Please read these instructions carefully before using the product.

UDI: Information on General safety and performance requirements is public available at <https://ec.europa.eu>
Basic UDI-DI: 4044508CeraEasy23


SSCP: The Summary of Safety and Performance is public available at <https://ec.europa.eu/tools/eudamed/#/screen/home>
If not in EUDAMED, the information shall be made available to the public upon request

Part A: Product Description
A1. Device Description

CERABLATE® easy/easy TC are a family of Legacy Devices (devices previously CE marked under directives 93/42/EEC, see MDCG 2020-6). They are 7F uni-directionally steerable quadripolar Ablation catheters for intracardiac applications. They are used for treatment of cardiac tachyarrhythmias, typically on the right side of the heart, by means of radiofrequency (RF) ablation therapy, for sensing of intracardiac signals and for diagnostic pacing.

The catheter can be connected to an ablation generator, as well as via the ablation generator with an electrophysiological measuring station.

For temperature sensing, a thermistor or thermocouple (TC) is integrated in the tip electrode (ablation electrode) and facilitates permanent temperature control, especially during RD ablation. The electrodes are made of platinum-iridium, a state-of-the-art material that has shown to have excellent conduction and thermal properties in the past 30 years of RF ablation. The catheters feature steerability of the distal tip together with torsion control that allow position of the catheter at the site of interest, which is typically the right atrium and/or ventricle. The deflection of the catheter tip via an internal pull-wire can be adjusted by operating the manipulator at the distal handle. With its special locking device it is possible to keep a selected curve in the desired curvature. It is possible to adjust the tightness of the locking device by each individual user.


Handle

The handle is designed with flattened areas similar to a screw-driver with a manipulator for curve adjustment. Below the usage and the performance characteristics are described.


Manipulator

The curvature of the catheter tip can be adjusted by moving the manipulator at the handle of the CERABLATE® easy/easy TC catheter forwards or backwards. The manipulator furthermore allows to increase or decrease the manipulator forces by turning the manipulator anti-clockwise or clockwise, respectively, and thus to regulate the fixation of the curvature.

Distal curvature control

Infinitely adjustable catheter curvature is achieved by moving the front part of the plastic handle (manipulator) forwards and backwards.

Adjusting the manipulation forces

Hold the handle tight and turn the manipulator according to the label or marking. Turning in the direction of the open lock icon reduces the manipulation forces. Turning in the direction of the closed lock icon increases the manipulation forces and latches the curve of the tip. The manipulator should be rotated by no more than 1 turn in each direction from the factory setting.

A2. Variants

There are different two curve shapes available Classic and Sigma in Diameters of 30/45/60 mm, to give physician a choice of selection. Smaller curves may be used for smaller hearts, however this is up to physician discretion.

Models:

CERABLATE easy® / easy TC®

Standard RF catheter

Curve form: Classic

Curve form: Sigma

Temperature sensor: Thermistor / thermoelement (thermocouple)

Classic

Sigma

Catheter shaft

The catheter shaft consists of a braided tube to facilitate torque transmission and kink resistance.

Distal part of the catheter

The distal part consists of individual poles enables precise recording of the intracardiac signals during mapping and RF dissipation. The steerability, flexibility and the possibility to vary the curvature of the distal part facilitate positioning and alignment of the catheter to the heart's anatomy.

Essential use information

Essential user information	Specification	Accuracy
Catheter Shaft size diameter	7F	7F (= 2.33 mm + 0.15 mm)

Catheter usable length	80 (Non-TC) 110 cm (TC+non-TC)	80 cm (-1/+5 cm) 110 cm (-1/+5 cm)
TIP size length	4 mm	4 mm ± 0.5 mm
TIP material	Platinum-Iridium	Platinum-Iridium
Ring electrode size length	2 mm	2 mm ±0.1 mm
Spacing	2-2-2 mm (TC+non-TC) 1-4-1 mm (non-TC) 2-5-2 mm (TC)	Tolerance: ±0.5 mm
Curves	30, 45, 60 mm	30 mm ± 0,3 mm 45 mm ± 0,5 mm 60 mm ± 0,6 mm

A3. Accessories

None

A4. Combinations with other medical devices

The following other medical devices are available separately:

Cat. No.	Name	For Catheter	Connection to Generator
82040	TX-HAT500-THR	Cerablate easy	OSYPKA HAT500 (THR)
82041	TX-HAT500	Cerablate easyTC	OSYPKA HAT500 (TC)
82042	TX-IBI	Cerablate easyTC	Abbott IBI T11/Ampere
82043	TX-SMARTABLATE	Cerablate easyTC	Stockert SMARTABLATE
82044	TX-QUIBIC/EP SHUTTLE	Cerablate easyTC	Stocker EP Shuttle/ Biotronik Qubic
82045	TX-HAT300/Atakr	Cerablate easyTC	OSYPKA HAT300/Atakr
82046	TX-nGen	Cerablate easyTC	Biosense Webster nGen

All RF ablation catheters are connected to an extension cable, which itself is then connected to an RF generator, this is state of the art. Through the extension cable, CERABLATE® easy/easy TC catheters are used in combination with other medical devices listed below:

- Radio-frequency ablation generators, e.g. OSYPKA HAT 300 Smart/HAT 500 or similar devices
- Via RF generator: External Pulse generators / Stimulators, e.g. Biotronik UHS 3000 or similar, ECG-monitors, e.g. Biotronik ECS 3000 or similar

A5. Intended Purpose

The Osypka CERABLATE® easy/Cerablate easy TC are a product family of temporary, uni-directional steerable multipolar catheters for intracardiac RF ablation as well as for sensing of intracardiac signals and diagnostic pacing. They are intended for the treatment of right-sided tachycardia termed atrio-ventricular nodal-re-entry tachycardia (AVNRT).

A6. Patient Population

The OSYPKA CERABLATE® easy/easy TC catheter for use in patients indicated for RF catheter ablation of right-sided cardiac tachycardia termed atrio-ventricular nodal-re-entry tachycardia (AVNRT).

A7. Intended Medical Indication (= Indication)

The Osypka CERABLATE® easy/Cerablate easy TC are intended for use in patients indicated for radio-frequency (RF) catheter ablation of right-sided cardiac arrhythmias termed atrio-ventricular nodal re-entry tachycardias (AVNRT). RF catheter ablation of cardiac arrhythmias, including AVNRTs is described in the corresponding guidelines and consensus statements from major cardiology societies, including the European Society of Cardiology.

A8. Expected clinical benefit

Allows termination of cardiac arrhythmias to enable normal life in sinus rhythm for patients.

A9. Contraindications

- Patients after acute myocardial infarction
- Systemic infections
- Acute local infections at puncture site
- Patient had shown sensitivity to foreign objects in the past or has an extreme allergy
- Insufficient or unsuitable tissue, e.g. tissue damaged by radiation or suppuration
- Physiological or anatomical abnormalities which could lead to post-operative complications, such as bleeding diathesis, reduced resistance against infections, etc.
- Hemodynamic instability
- Patients unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation
- Patients with prosthetic valves

A10. Possible complications/side effects

- Infection / sepsis
- Embolisms
- Coagulum formation, in particular with the use of catheters in power-regulated mode
- Arterial spasms (including coronary arteries)
- Perforations
- Pericardial effusion
- Pericardial tamponade
- Endocarditis
- AV-block
- Abnormalities of cardiac wall or heart valves
- Thrombosis of arteries or veins
- Obstruction of the left coronary artery
- Severe cardiac arrhythmias
- Ventricular fibrillation
- Stroke
- Death

A11. Warnings/Precautions

- Please read these instructions carefully before using the product.
- Use of the product in sterile conditions.
- The product may be contaminated after use. For this reason please observe your hospital's regulations regarding disposal.
- The product is not autoclavable. Single use only. Do not re-sterilise used products. Re-use brings the risk of infection and the possibility of malfunction!
- If the sterile packaging is being damaged or unintentionally opened before use, the product may be contaminated and must not be used.
- The product may only be used in a medical treatment facility that has been specifically set up for the appropriate application and by trained staff.

- The catheter should only be applied by specialists trained in the method of intracardial radiofrequency catheter ablation/electrophysiology and practiced in the technique.
- Observe the operating instructions of the used radiofrequency generator/ablation generator.
- The procedure must take place under X-ray conditions with appropriate X-ray equipment.
- It is imperative to carefully check the catheter and its accessories in connection with the radiofrequency ablation system before using it.
- Pay particular attention to whether the catheter shaft and the lead tip have been damaged in any way. Replace faulty products.
- The functioning of the manipulator should be tested before beginning.

Caution: The catheters are equipped with adjustable manipulation forces. The manipulator should not be turned too far against the specifications given in the product description in the direction of the open lock, since, otherwise, the buckling protection may get shifted and thus become ineffective. Products affected in this manner must be replaced immediately.

- The product must not be modified in any way.
- Do not advance or retract the catheter when there is distinct resistance. The cause of the resistance must first be determined radioscopically and appropriate steps must be taken to prevent patient injuries.
- The product may only be used if its safe application can be guaranteed. It is the doctor's responsibility to select the medically appropriate procedure and method
- Trouble-free functionality of the catheter can only be guaranteed in combination with the appropriate connection cables and external devices offered from Osyпка.
- Before using individual OSYPKA products or sets with non-OSYPKA products, the user must ensure the compatibility of the individual products with the specific application
- Use approved medicinal products only and apply them according to their operating instructions
- Observe the respective operating instructions when connecting and using external devices.
- When connecting the catheter to the ablation device, ensure compatibility of the connections with the ablation device. The catheter should only be used if the temperature display of the ablation device displays a realistic temperature
- Observe standard safety regulations applying to all electrical appliances.
- The catheters represent a direct, low impedance current path to the heart. It is known that small electric currents are sufficient to make the heart flutter. Do not touch the connector with your bare hand and do not allow it to make contact with electrically conducting or damp surfaces.
- Keep away any static electricity from the catheter system. Special attention must also be given to sufficient and central grounding of the operating table and the electrical devices (e.g. X-rays) used.
- Never touch the catheter tip and the indifferent lead at the same time, especially not during output.
- A working external pacemaker and defibrillator should be available before and during the entire electrophysiological examination. During defibrillation, defibrillation discharges can result in current flows through the pacemaker/catheter circuit. To protect the patient, the catheter circuit must therefore be interrupted during defibrillation.

Caution: If devices operated with mains voltage are used near the patient.

- When connecting the catheter to the ablation device, ensure compatibility of the connections with the ablation device. The catheter should only be used if the temperature display of the ablation device displays a realistic temperature
- The function of active implantable medical products such as pacemakers and defibrillators can be impaired by radio-frequency catheter ablation.
- Radio-frequency catheter ablation causes electrical discharges in the catheter tip/tool tip, which may influence the function of external devices. Observe the operating instructions of the surrounding devices.
- The catheter's function can be influenced by magnetic fields.

Caution: Do not apply the catheter during MRI examinations

- The application of external devices emitting electromagnetic radiation can lead to interference.
- The function of cardiological implants can be impaired by catheter/lead placement. Such systems could further impair the function of the ablation catheter. The patient's medical history must be taken into consideration.
- All improper use before or during the application, such as excessive bending or kinking of the catheter, can lead to the destruction of the control mechanism.
- Unforeseen complications may occur during transseptal access.
- Maximum duration of catheter use may not exceed 24 hours.
- Do not clean the handle with liquids. This could lead to a short circuit.
- If RF emission is output-controlled, this may lead to an increased complication rate. Thrombus formation and/or carbonization of the tissue in particular may increase. The possibility of the impairment of catheter function caused by overheating cannot be excluded in this case.

Caution: When using a defibrillator, disconnect the ablation catheter from the connection cable or remove the catheter from the endocardium. Do not defibrillate via the ablation catheter!

Caution: The use of external devices may cause unforeseen risks. Clarify these with the device manufacturer.

- When connecting CERABULATE® catheters to other ablation systems, first check for compliance of PIN allocation. If compliance with the PIN allocation cannot be guaranteed, CERABULATE® may not be connected to other ablation systems.

Caution: The products should only be operated in temperature-controlled mode. Use the output limit and impedance disconnection functions.

Part B: General Notes and Precautions
B1. Intended Users (=Use and handling)

Caution: Use of the product in sterile conditions.

The product may only be used in a medical treatment facility that has been specifically set up for the appropriate application and by trained staff. The product may only be used if its safe application can be guaranteed. It is the doctor's responsibility to select the medically appropriate procedure and method. The operating instructions are intended purely as general information on handling the product. The product must not be modified in any way.

Warnings / general information / precautionary measures must always be observed even if a differing procedure is chosen for medical reasons! Failure to observe these can result in complications and malfunction!

B2. Sterilisation and storage

The products are sterilised with ethylene oxide. Provided the sterilised product is stored correctly, it may be used up until expiry date (use before date) indicated on the package. Store the sterile product in its original packaging in a cool, dry place away from light at between 10°C and 25°C.

REUSE PRECAUTION STATEMENT

Not autoclavable. Single use only. Do not re-sterilise used product. Re-use bears the risk of infection and the possibility of malfunction!

B3. Disposal

The product may be contaminated after use. After use, this product may be a potential biohazard. Handle and dispose of in accordance with medical practice and applicable local, state and federal laws and regulations. For this reason please observe your hospital's regulations regarding disposal.

Any pointed or sharp objects such as needles, hypodermic needles or scalpels that were part or used with the product must be disposed of safely and separately in order to prevent injury and infection.

B4. Liability limitation

OSYPKA AG warrants its products to be free from defects in workmanship and materials under normal use. This warranty is valid only before the expiration date indicated on the product labeling. The product shall be used only in accordance with the labeling and shall not be modified or subjected to misuse, abuse, accident or improper handling. Proper handling and approved methods of use of the product are described in the „Instructions for Use“ document included with each product. OSYPKA AG disclaims all responsibility and liability for the use of its products in a manner not authorized or approved.

OSYPKA AG's liability under this warranty is limited to the replacement of its products. The foregoing warranty disclaims and is in lieu of all other warranties of merchantability or fitness for a particular purpose. OSYPKA AG disclaims all liability for incidental or consequential damages, losses, or expenses arising directly or indirectly from the use of this product, except as expressly provided by specific law.

OSYPKA AG does not assume or authorize any other person to assume any other or additional liability for loss, damage or expense in connection with this product. For further information, please refer to the OSYPKA AG General Terms and Conditions, which is available from OSYPKA AG or on the back of an OSYPKA AG invoice

B5. Reporting of serious incidents

The user of this product is obliged to report any serious incident that has occurred in relation to this device to OSYPKA AG and the competent authority in its country.

Part C: DIRECTIONS FOR USE

C1. Inspection of the product prior to use

The product, packaging and any accessories must be inspected for damage and malfunction before use. The product must not be used if any damage or malfunction is encountered. If the instructions for use are no longer available you can request replacement instructions from OSYPKA AG. The instructions for use may only be disposed of after all products in the product unit have been used.

C2. Unpacking the product

The catheter must be removed from its packaging under sterile conditions and used in a sterile work area.

C3. Catheter Insertion (Delivery Procedure)

- Insert the catheter into the vascular system via normal accesses for catheter systems.
- Connect the catheter to the radio-frequency ablation device (see operating instructions for the used device) or to the ECG recording device via suitable connection cables.
- Insert the catheter into the endocardial area using the X-ray and ECG system
- For better positioning in the heart, you can bend the catheter tip with the use of the control mechanism and thus change the catheter curvature. Information on the control mechanism is contained in Chapter 1 “Device Description”.

Caution: Maximum duration of catheter use may not exceed 24 hours.

Cleaning: Should you need to clean the catheters for any reason, be sure to protect the handle from penetrating fluids. Do not use solvents to clean the catheter.

Caution: Liquid in the handle will cause an electrical short circuit and may cause the catheter to malfunction.

C4. Medical Intervention

- You can start RF dissipation once the lead tip displays stable tissue contact with the intended ablation site.
- RF dissipation can be repeated at the same or other sites. Should the impedance increase for any reason, remove the catheter from the vessel. Clean the lead tip before repeating RF dissipation

C5. Removal Procedure

- Before pulling out the catheter from the patient, it must be checked with the help of X-ray that the tip of the catheter is not bent!
- After use, proceed in reverse sequence.
- Now remove the catheter according to the medical standard

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