CE 2797

Conformité Européenne (European Conformity), MDR 2017/745
 Conformité Européenne (Europäische Konformität), MDR 2017/745
 Conformité Européenne (Conformité Européenne), MDR 2017/745
 Conformité Européenne (Conformität europea), MDR 2017/745
 Conformité Européenne (Europäische Konformität), MDR 2017/745
 Conformité Européenne (conformidad europea), MDR 2017/745
 Conformité Européenne (Europæisk overensstemmelse), MDR 2017/745
 Conformité Européenne (Eurooppalainen vaatimustenmukaisuus), MDR 2017/745
 Conformité Européenne (Ευρωπαϊκή συμμόρφωση), MDR 2017/745
 Conformité Européenne (Atbilstība Eiropas standartiem), MDR 2017/745
 Conformité Européenne (Europese conformiteit), MDR 2017/745
 Conformité Européenne (Zgodność europejska), MDR 2017/745
 Conformité Européenne (Conformidade Europeia), MDR 2017/745
 Conformité Européenne (i överensstämmelse med EU-direktiv), MDR 2017/745
 Conformité Européenne (zhoda s požiadavkami EÚ), MDR 2017/745
 Conformité Européenne (Evropska skladnost), MDR 2017/745
 Conformité Européenne (Avrupa Uygunluğu), MDR 2017/745
 Conformité Européenne (Európai Megfelelőség), MDR 2017/745
 Conformité Européenne (Conformitate europeană), MDR 2017/745
 Conformité Européenne (Europos atitiktis), MDR 2017/745
 Conformité Européenne (Evropska skladnost), MDR 2017/745
 Conformité Européenne (Euroopa vastavus), MDR 2017/745
 Conformité Européenne (Европейско съответствие), MDR 2017/745

UK
CA 0086

United Kingdom Conformity Assessed

UKRP

United Kingdom Responsible Person:
 Dot Medical Ltd. trading at Silk Point,
 Queens Avenue, Macclesfield, SK10 2BB, United Kingdom

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


TME®

Temporary Myocardial Electrodes

Caution: Please read these instructions carefully before using the product.

UDI: Information on general safety and performance requirements is publicly available at <https://ec.europa.eu>
 Basic UDI-DI: 4044508TMEBV

SSCP: The Summary of Safety and Performance is public available at <https://ec.europa.eu/tools/eudamed/#/screen/home>
 If not in EUDAMED, the information shall be made available to the public upon request

Part A: Product Description
A1. Device Description

OSYPKA TME® temporary myocardial electrode is a device consisting of flexible insulated electrical conductors with one end connected to an external pacemaker pulse generator and the other end applied to the heart. The device is used to transmit a pacing electrical stimulus from the pulse generator to the heart and/or to transmit the electrical signal of the heart to the pulse generator. The conductors or wires are typically 60cm or 220cm long and can be fixated to the heart either with the help of a suture (fixation V and Loop), or directly with the attached heart needle (fixation type Z, T and S). The heart needle is cut off after fixation and disposed of. The thoracic needle is used to lead the wires to the outside of the patient’s body to allow connection to the external pulse generator. The thoracic needle(s) is then also disposed of.

OSYPKA TMEs are supplied as unipolar, bipolar bifurcated, bipolar, quadripolar and double-bipolar version. There are two options for connection to an external cardiac pacemaker or cable: with 2 mm adapter or Confix. The pre-assembled Confix connectors can be connected directly to an extension lead or to the cardiac pacemaker after the thorax needle has been cut off.

The device is supplied sterile, non-pyrogenic and intended for single use only.

Common Name: Pacing lead, pacing wire, heart wire, pacing electrode
 Device Name/Trade Name: TME® Temporary Myocardial Electrode

A2. Variants
Overview of the general features of OSYPKA TME® pacing wire types:

	Polarity	Lengths	Fixation Type	Thorax needle shape	Heart needle diameter	Connection with
OSYPKA TME	Unipolar, bipolar, quadripolar	60cm-250cm	Zigzag Tines V-hook Loop Straight	Curved or straight	0.30-0.46mm	2mm Adapter, Confix pin

A3. Accessories
 The following are accessories according to Article 2 (2) EU/VO 2017/745:

	Included in which type of TME sterile package			
Accessory	TME unipolar blue litz wire	TME unipolar white litz wire	TME bipolar	TME quadripolar
2mm Adapter blue	1x	-	1x	2x
2mm Adapter white	-	1x	1x	2x

A4. Combinations with other medical devices
 The TME® pacing wires are developed to operate with OSYPKA PACE101H, PACE 203H, PACE300 or DefiPace™ with 2mm adapters, Confix Pin and connection cables, such as OSYPKA D2-SP and OSYPKA D2-PSP connection cables. In the case of TME® with 2mm adapters, the adapters are included in the TME® package.

Cables			
Part Number	Product Name	Description	Length
81973	D2-SP	TME Connection Cable	2.5 m
81986BL	D2P-SP	TME Connection Cable blue	2.5 m
81986WS	D2P-SP	TME Connection Cable white	2.5 m
81975-3 81975-5 81975-10	D2-SP Adapter	D2-SP Cable	3 cm 5 cm 10 cm

External Devices		
Part Number	Product Name	Description
61216	OSYPKA PACE101H	External Pacemaker single chamber
61225 61225-ST	OSYPKA PACE203H	External Pacemaker dual chamber
61227 61227-ST	OSYPKA PACE300	External pacemaker, triple chamber, CRT
61228	DefiPace	External pacemaker, triple chamber, bi-atrial pacing & CV

A5. Intended Purpose
 OSYPKA TME® Temporary Myocardial Electrodes are intended for temporary cardiac stimulation and/or sensing of the cardiac chambers after cardiac surgery, in combination with an external pacemaker with or without connection cable.

A6. Intended Medical Indication (= Indication)
 OSYPKA TMEs are intended for temporary cardiac pacing and sensing after cardiac surgery to prevent postoperative bradycardia and conduction disorders.

A7. Patient Population
 OSYPKA TME® Temporary Myocardial Electrodes are intended for adults and children considered suitable by the physician. The use of OSYPKA TMEs is contraindicated when a permanent cardiac

pacemaker or internal cardioverter device (ICD) has already been implanted.

A8. Expected clinical benefit
 The OSYPKA TME® has the following clinical benefits:

- The measure of pacing/sensing of the post-operative patient allows the prevention of bradycardia and its deleterious side effects (e.g. hemodynamic instability) in the post-operative patient.

A9. Contraindications
 The use of OSYPKA TMEs is contraindicated when a permanent cardiac pacemaker or internal cardioverter device (ICD) has already been implanted. The particular medical condition and anatomy of the patient, however, may dictate the electrodes and implantation procedure to be used by the physician. OSYPKA TMEs may not be used for ventricular defibrillation or cardioversion.

A10. Possible complications/side effects

- Cardiac arrhythmias
- Bradycardia
- Perforation
- Pneumothorax
- Hemothorax
- Skeletal muscle and nerve stimulation
- Major bleeding / hemorrhage and tamponade especially during removal
- Infection, inflammation

A11. Warnings/Precautions

- OSYPKA TMEs may only be implanted by experienced, trained medical staff.
- TME® are ethylene-oxide sterilized prior to shipment. TME® must not be re-sterilized. Expired products shall be discarded.
- The packaging should be inspected prior to use. Should the seal or package be damaged, contact the local OSYPKA representative.
- Use of the product in sterile conditions.
- Pacing wires create a direct, low-resistance current path to the patient’s heart, possibly causing arrhythmias such as atrial fibrillation. Therefore, for devices that are used in the direct vicinity of the patient that operate on mains electricity it is vital to act in accordance with the instructions for electro-medical devices in order to prevent current leakage to the heart.
- Line-powered equipment used in the vicinity of the patient must be properly grounded.
- Avoid touching non-insulated parts of the product (e.g. in the process of connecting the non-insulated portion of the litz wire with the adapter) with bare hands or allow them to come into contact with electrically conductive or moist surfaces.
- Static electricity sources must be kept away from the stimulation system.
- Care must be taken to ensure that the operating table as well as any electrical equipment used is sufficiently and centrally earthed.
- Continuous, separate ECG monitoring of the patient is recommended when stimulating via OSYPKA TMEs since a significant increase in the stimulus threshold and a loss of sensation are possible.
- Equipment connected to TMEs shall comply with IEC 60601-1-1, type CF. Please follow the relevant instructions for use when connecting and using external devices (e.g. OSYPKA Pace203H, OSYPKA Pace300).
- The product is not suitable for ventricular defibrillation or electrical cardioversion. The use of this product for defibrillation with an external defibrillator may result in loss of function and scar formation.
- During defibrillation with an external defibrillator, OSYPKA TMEs must not be touched.
- Avoid damaging OSYPKA TMEs through the use of RF surgical equipment.
- OSYPKA TME® requiring 2mm adapters (included in the package) must not be connected to an external cardiac pacemaker without adapter under any circumstance. Stimulation or sensing loss may occur if only the exposed strand wire is connected.
- The needles must be handled with care and must not be bent any further.
- It is important that the electrodes are not placed poorly conducting fatty cardiac tissue or in tissue that is damaged, e. g. due to previous myocardial infarction.
- Patients with OSYPKA TMEs may not be exposed to significant electrical or magnetic fields (e.g. MRI). Exposure to magnetic fields may cause product and tissue heating, which may lead to failure of the product’s pacing and sensing function. Please refer to the 2021 ESC Guidelines on cardiac pacing for further information regarding handling of patients with pacing wires requiring an MRI : 2021 ESC Guidelines on cardiac pacing and cardiac resynchronization therapy (escardio.org). In the Guidelines a reference to OSYPKA TME is made through the publication by Homsí et al: Safety and Feasibility of Magentic Resonance imaging of the Brain at 1,5 Tesla in Patients with Temporary Transmyocardial Pacing leads. Thoracic and Cardiovascular Surgeon Vo, 67, No. 2/2019, pp. 86-91.
- Care should be taken such that the wires are not accidentally pulled on during transport of the patient.
- OSYPKA TMEs may remain in the body for a maximum of 10 days only.
- Leaving the Temporary Myocardial Electrodes implanted for longer than 10 days may result in difficulty or inability to extract the leads and/or bleeding.
- Removal of OSYPKA TMEs must be carried out with great care and should only be performed by an experienced physician.
- It is recommended that defibrillation equipment is readily available during temporary pacing lead insertion and removal.
- After removal patients should be monitored for at least 24 hours.
- TME® must not be cut off at the thorax, but should be completely removed from the patient. Product cut off at the chest may be prone to migration, infection or inflammation.
- The product is intended for single-use only. Possible harm due to infection/contamination and related consequences may occur to staff and/or patient if the product is re-used.
- Parts of the product may possibly contain cobalt residues. Possible carcinogenic side effects have to be considered, particularly for vulnerable patient groups such as children, pregnant and breast-feeding women. For more information please refer to: Kovovich M, Monnot A, Kougia D, More S, Wiley J et al. Carcinogenic hazard assessment of cobalt-containing alloys in medical devices: Review of in vivo studies, Elsevier Regulatory Toxicology and Pharmacology 122 (2021) 104910 Eichenbaum G, Wiley J, Fessel G, Qiu Q, Perkins L, et al, An integrated benefit-risk assessment of cobalt-containing alloys used in medical devices: Implications for regulatory requirement in the European Union, Elsevier Regulatory Toxicology and Pharmacology 125 (2021) 105004

Part B: General Notes and Precautions
B1. Intended Users (Use and handling)
 The product may only be used in a medical treatment facility that has been specifically set up for cardiac surgery and by a physician and surgical assistant specialized in cardiac surgery. The product may only be used if its safe application can be guaranteed. It is the doctor’s responsibility to select the medically appropriate procedure and method. The operating instructions are intended purely as general information on handling the product. The product must not be modified in any way. Warnings/ general information/ precautionary measures must always be observed even if a differing procedure is chosen for medical reasons! Failure to observe these can result in complications and malfunction!

B2. Sterilisation and storage
 The products have been sterilised with ethylene oxide. Provided the sterilised product is stored

correctly, it may be used up until expiry date (use before date) indicated on the package. Store the sterile product in its original packaging in a cool, dry place away from light at between 10°C and 25°C.

B3. Disposal

The product may be contaminated after use. After use, this product may be a potential biohazard. Handle and dispose of in accordance with medical practice and applicable local, state and federal laws and regulations. For this reason please observe your hospital's regulations regarding disposal. Any pointed or sharp objects such as needles, hypodermic needles or scalpels that were part or used with the product must be disposed of safely and separately in order to prevent injury and infection.

B4. OSYPKA AG Limited Warranty

OSYPKA AG warrants its products to be free from defects in workmanship and materials under normal use. This warranty is valid only before the expiration date indicated on the product labeling. The product shall be used only in accordance with the labeling and shall not be modified or subjected to misuse, abuse, accident or improper handling. Proper handling and approved methods of use of the product are described in the „Instructions for Use“ document included with each product. OSYPKA AG disclaims all responsibility and liability for the use of its products in a manner not authorized or approved.

OSYPKA AG's liability under this warranty is limited to the replacement of its products. The foregoing warranty disclaims and is in lieu of all other warranties of merchantability or fitness for a particular purpose. OSYPKA AG disclaims all liability for incidental or consequential damages, losses, or expenses arising directly or indirectly from the use of this product, except as expressly provided by specific law.

OSYPKA AG does not assume or authorize any other person to assume any other or additional liability for loss, damage or expense in connection with this product. For further information, please refer to the OSYPKA AG General Terms and Conditions, which is available from OSYPKA AG or on the back of an OSYPKA AG invoice.

B5. Reporting of serious incidents

The user of this product is obliged to report any serious incident that has occurred in relation to this device to OSYPKA AG and the competent authority in its country.

Part C: DIRECTIONS FOR USE

C1. Inspection of the product prior to use

The product, packaging and any accessories must be inspected for damage and malfunction before use. The product must not be used if any damage or malfunction is encountered. If the instructions for use are no longer available you can request replacement instructions from OSYPKA AG. The instructions for use may only be disposed of after all products in the product unit have been used.

C2. Unpacking the product

The product must be removed from its packaging under sterile conditions and used in a sterile working environment. After the outer sterile bag has been opened, the inner sterile bag must be carefully transferred to the sterile work area. Remove the TME® and the jig from the inner bag, detach the cardiac needle(s) (if applicable) from the jig, unwind the TME® carefully from the jig and detach the thorax needle.

C3. Positioning and attaching OSYPKA TME®.Z models

Fasten the distal end of the OSYPKA TME®.Z using the smaller, curved cardiac needle at the desired location of the heart in the myocardium (atrium, ventricle). Fastening is largely achieved as a result of the zigzag arrangement of the OSYPKA TME®.Z (Figs. 1a, 1b).

Do not use an additional ligature to fasten the zigzag anchoring attachment. Position the bipolar version in such a way that the hull (if applicable) is on and not in the myocardium (Fig. 1b). Cut off the distal end of the OSYPKA TME®.Z including the cardiac needle.

Fig. 1a: TME..Z, unipolar or bifurcated


Fig. 1b: TME..Z, bipolar


Fig. 1c: TME..Z in heart

C4. Positioning and attaching OSYPKA TME®.T models

Fasten the distal end of the OSYPKA TME®.T using the smaller, curved cardiac needle at the desired location of the heart in the myocardium (atrium, ventricle).

Fastening is largely achieved as a result of the anchors splaying and catching (Figs. 2a, 2b).

Splaying of the anchors is achieved by drawing the OSYPKA TME®.T through the cardiac muscle until the anchor fully emerges from the myocardium.

Subsequently carefully draw the OSYPKA TME®.T back until the distal sections of the anchors just protrude from the myocardium. The anchors will splay open in this process. Then draw the OSYPKA TME®.T again forwards slightly so that the anchors do not exert any pressure on the myocardium.

Do not use an additional ligature to fasten the anchoring attachment.

Position the bipolar version in such a way that the hull is on and not in the myocardium (Fig. 2b). Cut the distal end off leaving a clearance of at least 10 mm to the anchor.

(The 10 mm residual section and a narrow connecting bridge between the anchor and the proximal end ensure that the anchor does not shear off and remain in the heart when the OSYPKA TME®.T is removed.)


Fig. 2a: TME..T, unipolar or bifurcated


Fig. 2b: TME..T, bipolar


Fig. 2c: TME..T, in heart

C5. Positioning and attaching OSYPKA TME®.V models

Catch the V-end of the OSYPKA TME®.V in the purse-string suture. You may also set a U-stitch suture at any desired location in the myocardium and fasten the V-attachment of the OSYPKA TME®.V to it (see examples in Figs. 3a, 3b).


Fig. 3a: TME..V in heart

Fig. 3b: TME..V

C6. Positioning and attaching OSYPKA TME®.L (Loop) models

The OSYPKA TME®.L (Loop) is fastened by means of a ligature with a surgical thread (Fig. 4a). The suture should not be set too tightly in order to allow the OSYPKA TME®.L to be removed with a gentle pull after use. The OSYPKA TME®.L is removed as shown in Fig. 4b.

Caution: Do not draw the suture through the loop in order to avoid injury during removal (Fig. 4c).


Fig. 4a: TME..L in heart


Fig. 4c: TME..L

Fig. 4b: TME..L Removal

C7. Positioning and attaching OSYPKA TME®.S models

Fasten the distal end of the OSYPKA TME®.S using the smaller, curved cardiac needle at the desired location of the heart in the myocardium (atrium, ventricle).

Note: Depending on the surgeon's personal preference, the OSYPKA TME®.S can be anchored on or in the cardiac muscle in various ways with or without suture material (e.g. the exposed strand wires can be bent slightly with a scalpel after the cardiac needle is cut off).

C8. Drawing through the thorax

Insert the thorax needle through the thorax and draw the TME® out until no great tension occurs and there are no hindrances to fastening. The OSYPKA TME® can be fastened once more on the surface of the thorax.

Note: OSYPKA TMEs must be strain-relieved in order to ensure safe operation.

C9. Handling OSYPKA TMEs with Confix

Cut off the proximal residual section (thorax needle and wire) of the OSYPKA TME® direct after the proximal end of the Confix (see Fig. 5).

Confix connectors (= Confix) can be connected directly to an OSYPKA extension lead or to the OSYPKA external cardiac pacemaker.


Fig. 5: Cut off residual wire at Confix pin end

If the TME® is not connected to an external cardiac pacemaker or extension lead immediately, slip the protective caps (Confix holders) over the Confix connector (Fig. 6).


Fig. 6 Protective cap for Confix Pins

This protects the Confix from being touched.

C10. Handling OSYPKA TME® adapters (2 mm adapters)

TME® 2mm-adapters are to be used for all OSYPKA TME® models without Confix connectors (pins). Use one adapter for each strand wire: 1 adapter for unipolar OSYPKA TME® models, 2 adapters for bipolar OSYPKA TME® models and 4 adapters for quadripolar OSYPKA TME® models. Appropriate number of adapters are included in each package.

Cut off the thorax needle by cutting through the exposed litz wired at a distance of 1.5 cm from the insulation.

Caution: Make sure you only touch the OSYPKA TME® and the adapter pin with gloves.

Insert the individual litz wire strand through the hole of the adapter pin (Fig. 7a). Push the pin with the litz wire into the adapter (Fig. 7b). Pull the cap over the back end of the adapter (Fig. 7c). Now connect the adapter with an OSYPKA external pacemaker or cable according to the instructions for use. If the TME® is not in use, pull the adapter cap over the pin (Fig. 7d).

Caution: OSYPKA TME® adapters are only supplied, and can only be used with OSYPKA TME® electrodes without Confix.


Fig. 7a-7d: Attaching 2mm Adapters to the litz wire

C11. Connecting the TME® to an external cardiac pacemaker/extension cable

Connect the proximal end of the OSYPKA TME® featuring either a 2 mm Adapter or Confix pin directly, or with the use of an OSYPKA extension lead (e.g. D2-SP, D2P-SP), to an external OSYPKA cardiac pacemaker. Please follow the instructions for use of the external cardiac pacemaker and of the extension lead.

Caution: OSYPKA TMEs may only be used with external OSYPKA cardiac pacemakers (such as PACE 101H, PACE 203H, PACE 300) and the corresponding OSYPKA connecting leads. The instructions for use of the products concerned must in all cases be followed.

Caution: OSYPKA TMEs may not be connected to the external cardiac pacemaker with exposed strand wires – without adapter or without Confix – under any circumstance. (Stimulation or sensation loss may occur if only the exposed strand wire is connected.)

Note: After connecting the OSYPKA TME® and when the thorax is still open, please check the position of the attachment, the sensation and stimulus threshold. If it is not possible to successfully generate stimulation or sensation thresholds, connect the OSYPKA TME® with the reverse polarity to the external cardiac pacemaker. Use a new OSYPKA TME® at a new location on the heart if sensing/ stimulation is subsequently still not possible.

C12. Post-operative pacing

Postoperative pacing should be facilitated according to the current guidelines and according to the IFU of the external devices.

Note: Pacing and sensing thresholds should be checked daily and adjusted, if necessary.

C13. MRI examinations

Note: The product must NOT be used during MRI examinations. Exposure to magnetic fields may cause product and tissue heating, which may lead to failure of the product's pacing and sensing function. Please refer to the 2021 ESC Guidelines on cardiac pacing for further information regarding handling of patients with pacing wires requiring an MRI : 2021 ESC Guidelines on cardiac pacing and cardiac resynchronization therapy (escardio.org).

C14. Care when not in use

If the patient does not require any more cardiac stimulation for a temporary period, the OSYPKA TMEs can be disconnected from the external cardiac pacemaker. When not in use, the TMEs can easily be fastened at the patient's skin with a medical strip to avoid tangling of the litz wires. The connections must be protected from contact by attaching the protective caps.

Note: In TME® models with 2mm adapter: Do not remove the 2mm Adapters from the wire!

Caution: Make sure you only touch the adapters or Confix of the OSYPKA TMEs with gloves. The electrostatic potential between the user and the patient should be equalized before handling the external cardiac pacemaker or the OSYPKA TMEs, for example by touching the patient at some distance from the OSYPKA TMEs.

C15. Removing OSYPKA TMEs

When the patient no longer requires any further cardiac stimulation, remove the OSYPKA TMEs by pulling carefully on the ends of the TME®. TMEs must not be cut off at the thorax. After removal, patients should be monitored for at least 24 hours.

Caution: OSYPKA TMEs must be removed after 10 days at the latest. Complications may arise if the TMEs are remained in situ for more than 10 days.

Caution: Please note the following approach to removing multipolar OSYPKA TME® models in order to avoid irritation during removal: Gently pull each individual litz wire a few centimetres out in order to release the attachment to the cardiac muscle. The OSYPKA TME® can then be drawn out completely.

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